



RoWalker


IMPROVE PATIENT MOBILITY AND PREVENT FALLS

RoWalker

MOBILIZE PATIENTS EARLY AND SAFELY



- 400 lb weight capacity
- Mobilize Patients Early and Reduce Length of Stay
- Challenge patients to walk farther
- The RoWalker carries all the equipment needed for complex patients
- Reduce risk of falls



IV POLE

HANDGRIPS

PADDED
ARMRESTS

HEIGHT
ADJUSTABLE

ACCESSORY
BASKET

SAFETY SEATS
- Flip up and
Rotate Out-

OXYGEN
HOLDER

DIRECTIONAL
WHEEL LOCK

BRAKES



Twist knob below platform to widen or narrow platform width as desired



Release the hydraulic latch to allow platform adjustment. Close latch once platform is in desired position



Apply directional locks on both sides to keep device moving in straight path, apply directional lock to one side to correct a gait drift.

Contact SPH Medical

Tel: 1-844-377-4633 | Email: info@sphmedical.com | Web: www.sphmedical.com